




- Ⓓ **Montageanleitung**
- ⒼⒷ **Mounting instructions**
- Ⓕ **Notice de montage**
- ⓃⓁ **Montage-instructie**
- ⓲ **Istruzioni di montaggio**
- Ⓔ **Instrucciones de montaje**
- Ⓟ **Instrukcja montażu**

D	Inhalt			I	Contenuto		
	Wichtige Hinweise	3			Informazioni importanti	7	
	Gebrauchsanleitung	10			Istruzioni per l'uso	10	
GB	Content			E	Contenido		
	Important information	4			Indicaciones importantes	8	
	Instructions for use	10			Instrucciones de uso	10	
F	Sommaire			PL	Treść		
	Remarques importantes	5			Ważne wskazówki	9	
	Notice d'utilisation	10			Instrukcją obsługi	10	
NL	Inhoud						
	Belangrijke tips	6					
	Gebruiksaanwijzing	10					


950.51.10090
950.51.10190

950.51.11090
950.51.11190

D

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


Purpose

The HEWI suspended seat is a seat for people with and without physical disabilities. The HEWI suspended seat is designed for a maximum vertical static load of 150 kg.

Mounting information

The following points must be observed:

- The suspended seat is suspended from the appropriate HEWI handrail (dia. 33 mm) from above. The support separation must not exceed max. 850 mm (horizontal bar).
- Ensure that the wall structure and the handrail fittings are able to permanently withstand the anticipated load. Lightweight wall constructions must have a suitable supporting framework.

Warning and hazard notices

- Improper mounting or use can lead to accidents involving serious injuries.
- Do not misuse the suspended seat, e.g. as a climbing aid.
- A certain vertical incline is normal when the seat is under load. In the event of damage, the suspended seat may only continue to be used following inspection by a specialist company.

Care and cleaning instructions



HEWI products are easy-care. Generally speaking, an occasional wipe with a damp cloth is enough. When using detergents, please observe the following:


- Only use detergents with a pH of between 6 and 8.
- Detergents containing active substances such as phosphates, soap and tensides may be used.
- Follow the detergent manufacturer's instructions (concentration, contact time).
- Detergents containing acids, alkalis, bleach or scouring agents may not be used.
- Do not use scouring cleaning equipment, such as brushes

Exclusion of liability

HEWI accepts no liability whatsoever for damage to persons or property that occurs or has occurred as a result of improper mounting or maintenance, particularly if the safety instructions and other points to note listed are not or have not been observed in full. HEWI accepts no liability whatsoever for damage occurring from incorrect cleaning.

For further information, please contact HEWI Technical Support, tel. +49 5691 82-300.

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HEWI Heinrich Wilke GmbH